You are not alone!

Join Us!

A community of women
for support and advice.

    @IntimateRose

Dr. Amanda Olson, DPT, PRC

VELVET ROSE
Personal Lubricant

The Essentials

Perfect with the Dilators

ENCHANTED ROSE
Organic Feminine Balm

Shop @:
IntimateRose.com

Manufactured For: Intimate Rose
North Kansas City, MO 64116
Made in PRC

Rectal Dilators

Includes English, French, German, Italian, Polish, Portuguese,
Spanish, Swedish, Dutch, Turkish, and Czech Translations.

Keys to Success

The keys to success are consistency and routine practice. Daily use of trainers coupled with relaxation techniques will result in achieving your goals.

How to Start

1. Wash the product with warm water and soap.
2. Designate a place in your home that is comfortable for your training. Choose a place that is quiet and comfortable.
3. Use a generous amount of water-based lubricant on the product and yourself. The use of a water-based lubricant is important to preserve the medical grade silicone of the trainer.
4. Select a size that looks appropriate for you.

5. Start by lying on your back with your knees bent and feet planted. Some people may prefer to lie on their side instead; if that's the case, be sure to have your knees bent and the top leg supported by a folded pillow between your knees.

6. Begin your training by breathing in and allowing your belly to expand, followed by exhaling, allowing the belly to slowly fall. Gently and carefully place the first ½ inch to 1 inch while maintaining a steady breathing pattern. Then gently place the rest of the trainer.

7. Repeat the slow, deep breathing cycle for 15 minutes.

8. Repeat this process 1-2 times per day, and progress to larger sizes until you can tolerate your goal size without pain.

General Care and Maintenance

It is important to clean with warm water and soap before and after each use. Store in a cool, dry place and do not share them with other people. If you find a defect, contact us for a replacement at Support@IntimateRose.com.

To maintain your progress after reaching your goals, it may be helpful to continue your 15-minute breathing practice up to 4-6 times per week for several weeks to months. You may find that once you have trained your body, you no longer need to practice!

Warning:

Non-sterile device is designed for single person use only. Latex free. Do not use silicone-based lubricants. Use slow, gentle movements.

Indications:

The trainers may be used to prepare for an examination, medical procedure, or relief from stenosis.

Contraindications:

Acute genital tract infections and pelvic infections. Ask your doctor before using when pregnant.

Important Notes Regarding Bleeding:

Spotting after surgery or radiation treatment is normal. If heavy bleeding occurs contact your physician.

Important notes on hygiene:

Wash hands and trainers with mild soap and warm water prior to use. Rinse thoroughly. If trainer feels cold, hold it under warm tap water and dry it well. Apply a generous amount of Intimate Rose Velvet Rose Lubricant.

Recommendations:

1. It is recommended that you see a pelvic floor physical therapist or occupational therapist for one-on-one guidance and a plan of care to help you achieve your goals.

2. Pelvic floor stretches are helpful in conjunction with training for improved pelvic floor function. Videos and guides on these exercises can be accessed for free on **IntimateRose.com** under "Guides."

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