You are not alone!
Join Us!

A community of women
for support and advice.

Dr. Amanda Olson, DPT, PRPC

VELVET ROSE
Personal Lubricant

The Essentials
Perfect with the Pelvic Wand

ENCHANTED ROSE
Organic Feminine Balm

Shop @:
IntimateRose.com

Manufactured For: Intimate Rose
North Kansas City, MO 64116
Made in PRC

**Bendable
Pelvic Wand**

Includes English, French, German, Italian, Polish, Portuguese,
Spanish, Swedish, Dutch, Czech and Turkish Translations.

Keys to Success

The Bendable Wand is designed to massage the pelvic floor through gentle pressure. Regular use of the wand, along with practicing relaxation of pelvic muscles, can aid in easing pelvic discomfort. For personalized advice, consult your healthcare provider.

How to Start

- 1. Wash Up:** Wash hands and the wand with mild soap and warm water before use.
- 2. Select End:** Both ends can be used in the vagina or rectum. The pointed end may be helpful in relieving tender points near the entrance of your body, while the thinner end may be more helpful in reaching the deeper pelvic muscles inside.
- 3. Apply Lubricant:** Use plenty of water-based lubricant on the wand and your body.
- 4. Get Comfortable:** Lie on your back with knees bent and feet planted; or lie on your side with a pillow between your knees.
- 5. Insert the Wand:** Begin by breathing in and allowing your belly to expand, followed by exhaling, allowing your belly to slowly fall. The act of slowly exhaling helps to naturally relax the pelvic floor muscles. Gently insert the wand.
- 6. Find Tender Points:** After insertion, press the wand into the internal tissue and move it clockwise until you find a tender spot. Press gently with the same amount of pressure you would use to check a tomato for ripeness. Don't squish your tomato!

7. Apply Pressure: Hold the wand on the tender spot and move your knee side to side until the pain eases. Hold for 1-2 minutes, then release. Continue to breathe!

8. Repeat: Continue to work clockwise searching for tender spots, and do the same exercise above each time one is found. Do this 1 or 2 times a day as needed.


Don't Forget the Lube!

Water-based lubricant is recommended for use with the pelvic wand. Explore Intimate Rose FDA-cleared, water-based options, original Velvet Rose and glycerin-free Serenity.

Bendable Wand Positions

Explore various wand configurations and body positions, such as lying on your back with knees bent, on your side, or with legs supported by a couch or chair, to find what works best for you.

Find More Detailed Tutorials at
intimaterose.com/bendable

Additional Resources

Vaginal and rectal dilators often complement the use of pelvic wands for specific conditions.

Whether for hormone balance, PCOS management, bladder support, or vaginal pH restoration, learn why 100,000+ women trust Intimate Rose Supplements!

Discover Enchanted Rose, a natural, daily-use vaginal balm to soothe and moisturize dry or irritated tissue.

Visit IntimateRose.com to learn more!

INDICATIONS:

Use to help massage the pelvic floor.

CONTRAINDICATIONS:

If pregnant, please check with your doctor. Do not use wand if you have active

vaginal or genital infections, undiagnosed bleeding, or are currently placed on pelvic rest*.

*Pelvic rest is common soon after having a vaginal surgery, cervical LEEP, cone procedure, or giving birth.

GENERAL CARE:

- Clean your hands and the wand with warm water and soap before and after each use.
- Store in a cool, dry place.
- Intimate Rose pelvic wands are not made with natural rubber latex.

WARNINGS:

- Prolonged use may cause vaginal abrasion or ulceration.
- If symptoms, such as bleeding, vaginal discharge, unpleasant odor, pain, or severe itching occur, discontinue use and consult your physician immediately.
- Do not share with others. Non-sterile device is designed for single person use only.
- Use water-based lubricants only. Do not use silicone-based lubricants.

PRECAUTIONS:

Before each use, check your device for signs of wear/damage (discoloration or splits/cracks/pitting). If signs of wear or damage are noted, discontinue use. Contact us for a replacement at

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