You are not alone!

Join Us!

A community of women
for support and advice.

@IntimateRose

Dr. Amanda Olson, DPT, PRPC

VELVET ROSE
Personal Lubricant

The Essentials

Perfect with the Dilators

ENCHANTED ROSE
Organic Feminine Balm

Shop @:
IntimateRose.com

Manufactured For: Intimate Rose
North Kansas City, MO 64116
Made in PRC

Dilators

**Includes French, German, Italian, Polish, Portuguese,
Spanish, Swedish, Dutch, and Turkish Translations.**

Keys to Success

To decrease pain with vaginal penetration, the keys to success are consistency and routine practice. Daily use of the vaginal trainers, coupled with relaxation techniques and focused attention on training the muscles, will result in achieving your goals. Your health care provider may have a unique training plan for you outside of the recommendations made here.

How to Start

1. Wash the dilators with warm water and soap.
2. Designate a place in your home that is comfortable for your dilator practice (also referred to as vaginal trainer). Choose a place that is quiet and comfortable.
3. Use a generous amount of water-based lubricant on the vaginal trainer and the opening of the vagina. The use of a water-based lubricant is important to preserve the medical grade silicone of the vaginal trainer.
4. Select a vaginal trainer that looks appropriate for you or that a health care provider has suggested.

5. Start by lying on your back with your knees bent and feet planted. Some people may prefer to lie on their side instead; if that's the case, be sure to have your knees bent and the top leg supported by a folded pillow between your knees.

6. Begin your training by breathing in and allowing your belly to expand, followed by exhaling, allowing the belly to slowly fall. Gently bring the vaginal trainer to the opening of the vagina and carefully insert it on an exhale.

7. Keep the vaginal trainer inserted in the vagina and repeat the slow, deep breathing cycle for 15 minutes or the length of time recommended by your health care provider.

8. If you experience discomfort, visual imagery is helpful for relaxing the pelvic floor muscles around the opening of the vagina. There are several images that can work. Pick the one that is best for you: "imagine that the vagina is like a rose, blooming outward and opening," "imagine that the vagina is like an umbrella, opening and expanding," or "imagine that the vagina is like an elastic band, flexible and mobile."

9. Repeat this process 1-2 times per day, and progress to larger trainers until you are able to tolerate a dilator of your goal size pain free.

General Care and Maintenance

It is important to clean the dilator with warm water and soap before and after each use. Store dilators in a cool, dry place and do not share them with other people. If you find a defect in your dilator, contact us for a replacement at Support@IntimateRose.com.

For general maintenance of your training, once you have reached your end goal, it may be helpful to continue your 15-minute breathing and vaginal trainer practice up to 3 times per week for 3-4 weeks. Your health care provider may also give you a specific maintenance plan. You may find that once you have trained your body, you no longer need to practice with your trainer!

Warning:

Non-sterile device is designed for single patient use only. Latex free. For use under the recommendation of a health care provider. Do not use oil or silicone-based lubricants. Use Intimate Rose's Velvet Rose Lubricant. Do not force the dilator into the vagina. Use slow, gentle movements.

Indications:

The dilators may be used to enlarge the vagina for an examination, surgical procedure, or difficulties with penetration for intercourse.

Contraindications:

Dilators are contraindicated in acute genital tract infections and pelvic infections. A physician must advise a pregnant patient on use of this product.

Important Notes Regarding Bleeding:

Spotting after surgery or radiation treatment is normal. If heavy bleeding occurs with use of the dilator, contact your physician.

Douching is not necessary.

Important notes on hygiene:

Wash hands and dilators with mild soap and warm water prior to use. Rinse thoroughly. If dilator feels cold, hold it under warm tap water and dry it well. Apply a generous amount of Intimate Rose Velvet Rose Lubricant to the end of the dilator and insert slowly to restore vaginal elasticity.

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